
Withings BPM Connect

INSTRUCTIONS MANUAL

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EN | FR | DE

THANK YOU FOR CHOOSING WITHINGS BPM CONNECT

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This guide is available in additional languages at: www.withings.com/guides

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
1

Visit go.withings.com from your smartphone or tablet, and follow the instructions to install the Health Mate app.

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2

In the app, follow the instructions to install Withings BPM Connect.

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How to take a blood pressure measurement

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
2

Unroll cuff. Place your arm inside it.

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3

Tighten the cuff around your arm. The tube should be positioned against your inner arm.

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
4

Place your arm on a table and level with your heart.

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10

5

Turn on your blood pressure monitor to start the measurement. You can select a simple or triple measurement with a long press on the button. Press the button again to start the measurement.

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
12

Requirements

In order to use your BPM Connect, you need an iOS (10.0 or higher) or Android (6.0 and higher) device to install it. Thereafter, the product can be used in stand-alone, thanks to the Wi-Fi and Bluetooth® connections.

Important safety information

Consult your doctor during pregnancy, or if you suffer from arrhythmia and arteriosclerosis. Please read this section carefully before using the Wi-Fi Smart Blood Pressure monitor.

Intended use

The BPM Connect is a digital monitor intended for use in measuring blood pressure, heart rate. The device is intended to be used in a human adult population with an arm circumference between 9 inches to 17 inches (22 cm to 42 cm).

The BPM Connect is a medical device.

Contact your physician if hypertensive values are indicated.

General safety and precautions

- Do not forcibly bend the arm cuff.
- Do not inflate the arm cuff when it is not wrapped around your arm.
- Do not apply strong shocks and vibrations to the blood pressure monitor or drop it.
- Do not take measurement after bathing, drinking alcohol, smoking, exercising or eating.
- Do not immerse the arm cuff in water.
- Do not use with a pacemaker, a defibrillator or other electric implant.
- Use on adults only.
- Do not use on children or pets.

Android is a trademark of Google Inc.

iOS is a trademark of Apple Inc.

Caution

General usage

Always consult your doctor. Self-diagnosis of measurement results and self treatment are dangerous. People with severe blood flow problems, or blood disorders, should consult a doctor before using the blood pressure monitor. Cuff inflation can cause internal bleeding. Operational factors such as common arrhythmias, ventricular premature beats, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, pre-eclampsia or renal disease can affect the performance of the automated sphygmomanometer and/or its blood pressure reading. This device is a precision measuring equipment liable to be understood by lay user but it still should be handled with care. A long exposure of the device to lint, dust or sunlight might reduce its life time or damage it. Damaged cuff or sensor might lead to incorrect measurements. No usage for USB port. The user is an intended operator. Parts in contact with the skin: cuff. Measurements can be affected by extreme temperatures, humidity & altitude.

- Do not leave the blood pressure monitor unattended with infants or persons who cannot express their consent.
- Do not use the blood pressure monitor for any purpose other than measuring blood pressure.
- Do not disassemble the blood pressure monitor.
- Do not operate the blood pressure monitor in a moving vehicle (car, airplane).
- Do not use a cellular phone near the device.
- Do not use the device with the USB cable plugged.
- Improper continuous pressure of cuff or too frequent measurements may interfere blood flow and result harmful injury to the users. Check if the use of the device does not result in prolonged impairment of your blood circulation.
- Do not apply the cuff over the users' arm bearing a wound or medical treatment and so on, as this can cause further injury.
- Cuff should be used on the arm on the side of a mastectomy.
- Use of the device can temporarily cause loss of function of simultaneously used equipment on the same limb.

Cleaning

- Do not use an alcoholic-base or solvent agent to clean the device
- Clean the device with a soft and dry cloth
- The dirt on the cuff can be cleaned by a moisten cloth and soap
- Do not flush the device and cuff with much water
- Do not dismantle the device or disconnect the cuff or try to repair by yourself. If any problem happens, refer to the distributor
- Do not operate the device under severe environment of extreme temperature or humidity, or direct sunshine
- Do not shake the unit violently
- Do not submerge the device or any of the components in water
- Do not let the device under strong shocks, such as dropping the unit on the floor

Storage

- Store the device and the components in a clean and safe location
- If storage conditions are different from the usage conditions indicated in this document, please wait 30 minutes before taking a measurement.

Maintenance

If you can't fix the problems using the troubleshooting instructions, request service from your dealer. Manufacturer will make available on requested circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist manufacturer's staff or authorized representative for repair. It is generally recommended to have the device inspected every 2 years, to ensure proper functioning and accuracy. Do not use the device while doing maintenance steps.

Error messages

15

Error message	Causes	Countermeasures
Measurement could not be performed. Please try again. If the problem occurs again, please contact customer service.	Hardware auto-check error.	Measurement could not be performed. Please try again. If the problem occurs again, please contact customer service.
Connection error. Quit the application and reconnect your blood pressure monitor.	Communication error between the Blood Pressure Monitor and the collateral device.	<ol style="list-style-type: none">1. Push the blood pressure monitor button for 3 seconds to turn it off. Quit the application and push the blood pressure monitor button to turn it on.2. If the problem occurs again, please contact customer service.
Please wait until the cuff is totally deflated before starting a new measurement. Stay still during measurement.	<ul style="list-style-type: none">- Cuff is not deflated before measurement.- Noise interference.- Excessive user movements.	<ol style="list-style-type: none">1. Please wait until the cuff is totally deflated before starting a new measurement. Stay still during measurement.2. If the problem occurs again, please contact customer service.

Error messages

16

Error message	Causes	Countermeasures
Check that the Blood Pressure Monitor is correctly positioned on your arm and that measurement is performed in good conditions.	<ul style="list-style-type: none">- User movement, or not relaxed, or talking ...- Cuff is not correctly fastened- Cuff broken- Pump or valve failure- Pressure overflow- Noise interference- User movement, or not relaxed, or talking ...- User special characteristics- Battery level is too low	<ol style="list-style-type: none">1. Check that the Blood Pressure Monitor is correctly positioned on your arm and that measurement is performed in good conditions.2. Low battery level. Load your BPM Connect with a micro USB cable.3. If user has special characteristics, please contact your physician.4. If the problem occurs again, please contact customer service.
Low battery level.	Battery level is too low.	Load the device with the provided charging cable.

Troubleshooting

17

Problem	Remedy
The inflation action cannot be done or the air pressure cannot rise	<ol style="list-style-type: none">1. Check the cuff position and fasten the cuff correctly and remeasure the pressure again2. Check the blood pressure monitor connection to the iOS or Android device
The low battery icon is displayed on the device	Reload the device with the cable provided
The blood pressure cannot be taken and the application shows an error message	<ol style="list-style-type: none">1. Re-fasten the cuff2. Relax yourself and sit down3. Keep the cuff and heart at the same level during the measurement period4. Keep silent and still during measurement5. If the patient has severe heart beat problem, then the blood pressure may not be read correctly
Under normal measuring circumstance, the reading at home is different from that of the clinics or each measurement has various reading	<ol style="list-style-type: none">1. The variation is due to the different environments2. The blood pressure is changing according to the physiological or psychological status of the human body3. Show your recorded values to your physician

Specifications Technical data

Product description:	Digital automatic Blood Pressure Monitor
Model:	WPM05
Blood pressure measurement method:	Cuff oscillometric method
Cuff inflation:	Automatic inflation with air pump at 15 mmHg/s
Pressure sensor:	Gauge sensor
Measurement range (pressure):	0 to 285 mmHg, DIA 40 to 130 mmHg, SYS 60 to 230 mmHg
Measurement range (pulse):	40 to 180 beats/min
Accuracy (pressure):	Within +- 3 mmHg or 2% of reading
Accuracy (pulse):	Within +-5% of reading
Sensor:	Semiconductor pressure sensor
Operating conditions:	5 to 40°C, 15 to 90% RH, atmospheric 86Kpa-106kpa, altitude: 2000m
Storage and transport conditions:	-25 to 70°C, 10 to 95% RH, atmospheric 86Kpa-106kpa, altitude: 2000m

Arm type:	Use on left arm.
Power source:	5V 1A
Weight:	Approx. 245g
Accessories:	Charging cable, Instruction manual
Product life:	3 years
Wireless transmission:	Wi-Fi and BLE

Note: Specifications are subject to change without prior notice or any obligation on the parts of the manufacturer

Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal.

The symbol applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste.

At the end of devices useful life, the user must deliver it to the able collecting centers for electric and electronic garbage, or give back to the retailer when purchasing a new device. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the environment and health. In case of abusive disposal of device by the user, will be applied administrative endorsements in compliance with current standard. The device and its parts is mared with regard to disposal, as appropriate, in accordance with national or regional regulations.

Version 1.0 April 2019

FCC Statement

Federal Communications Commission (FCC) Statement 15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) this device may not cause harmful interference and
- 2) this device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

RF Statement

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following. | Interference may occur in the vicinity of equipment marked with | Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. | The use of accessories and cables other than those specified may result in increased emissions or decreased immunity | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter | The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Any other accessories, transducers and cables may result in increased emissions or decreased immunity and EMC performance. | The device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, which should be observed to verify normal operation in the configuration in which it will be used. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following. Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

Guidance and manufacturer's declaration-electromagnetic emissions

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
CE emissions CISPR11	Group 1	The Wi-Fi Smart Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RE emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	The Wi-Fi Smart Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	

Declaration – electromagnetic emissions and immunity

23

for equipment and systems that are not life-supporting and are specified for use only in a shielded location

The Wi-Fi Smart Blood Pressure Monitor declaration electromagnetic immunity

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RFIEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A N/A	N/A
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol.

Declaration – electromagnetic emissions and immunity

for equipment and systems that are not life-supporting and are specified for use only in a shielded location

Declaration – electromagnetic immunity

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	N/A N/A	A Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	-5 % UT(95% dip in UT) for 0.5 cycle, -40 % UT (60 % dip in UT) for 5 cycles, -70 % UT(30 % dip in UT) for 25 cycles, -5 % UT (95 % dip in UT) for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Withings Two (2) Years commercial guarantee Withings BPM Connect | Wi-Fi Smart Blood Pressure Monitor

Withings SA, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, France ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

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This device complies with Part 15 of the FCC Rules and with Industry Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications made to this equipment not expressly approved by Withings may void the FCC authorization to operate this equipment.

Radiofrequency radiation exposure information: The radiated output power of the device is far below the FCC radio frequency exposure limits. Nevertheless, the device shall be used in such a manner that the potential for human contact during normal operation is minimized.

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56

This device complies with the following normative documents : COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 concerning medical devices as amended by Directive 2007/47/EC

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
EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes. Reference to standards contd. | EN ISO 14971: Medical devices - Application of risk management to medical devices. | IEC/EN 60601-1: Medical electrical equipment - Part 1: General requirements for basic safety and essential performance. | EN 1060-3: Non-invasive sphygmomanometers, Part 3: Supplementary requirements for electromechanical blood pressure measuring systems. | EN 1060-4: Non-invasive sphygmomanometers. Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers. | IEC/EN 60601-1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. | IEC 80601-2-30: Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers. | EN 300 328: Electromagnetic compatibility

and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the Directive (2014/53/EU). | EN 301 489-1: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements. | EN 301 489-17: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements | IEC/EN60601-1-2: Medical electrical equipment: Part 1-2: General requirements for basic safety and essential performance-collateral standard electromagnetic compatibility | EN 55011: Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement. | FCC part B 15B: Electromagnetic Compatibility. | FCC Rule Part: 15.247 Cat: DSS (Bluetooth) FCC Rule Part: 15.247 Cat: DTS (BT4.0). | EN ISO 10993-1: Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process. | EN ISO 10993-5: Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity | EN ISO 10993-10: Biological evaluation of medical devices. Tests for irritation and skin sensitization | IEC/EN 60601-1-6: Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability | ANSI/AAMI/ISO 81060-2 Non-invasive sphygmomanometers Part 2: Clinical validation of automated measurement type | EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

Manufacturer/Fabricant/Hersteller:

WITHINGS SA,
2 rue Maurice Hartmann,
92130 Issy-les-Moulineaux, FRANCE
+33141460460

IP22

Ingress of water or particulate matter

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Complies with waste electrical and electronic equipment directive

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Frequency band :
2402 - 2480 MHz
Maximum output power : 4.45 dBm

EN Withings SA hereby declares that the device Withings BPM Connect is in conformity with the essential requirements and other relevant requirements of Directive 2014/53/EU.

The declaration of conformity can be found at: withings.com/compliance.

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Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Need help?

Besoin d'aide? | Brauchen Sie Hilfe?

→ withings.com/support

IM_WPM05_ENFRDE_001

withings

Withings BPM Connect | Wi-Fi Smart Blood Pressure Monitor | V1.0 - Apr.2019

© 2019 Withings SA. All rights reserved.