
Portable Mesh Nebulizer NE-U300 Instruction Manual

3156785-2A
IM-NE-U300-EO-01-11/2025
Doc: LFS-EN-UM-17C-EO
Rev.: A/0
Issue date: 2025-11-11

EN

DE

FR

ES

IT

NL

RU

TR

PL

SV

AR

Table of Contents

1. Intended Use -----	4
2. Contraindications -----	5
3. Important Safety Notes -----	5
4. Working Principle -----	6
5. Package Contents and Overview -----	7
6. Product Technical Parameters -----	8
7. Installation and User Instructions -----	14
8. Cleaning and Disinfection -----	22
9. Storage and Maintenance -----	28
10. Trouble Shooting -----	29
11. Note & Warning -----	31
12. EMC declaration -----	37
13. Signs and Symbols -----	43
14. Warranty -----	46
15. Configuration List -----	47
16. Disclaimer Clause -----	48

- Thank you for purchasing the OMRON Portable Mesh Nebulizer.
- Be sure to read this instruction manual carefully before using this device, so that you can use it safely and correctly.
- Please keep this instruction manual in a place convenient for reading at any time.
- The illustrations in this instruction manual are schematic diagrams.

1. Intended Use

1.1. Intended purpose

The Portable Mesh Nebulizer is intended to convert certain inhalable drugs into an aerosol form for inhalation by a patient for the treatment of obstructive airway diseases (OADs), such as asthma, chronic obstructive pulmonary disease (COPD), and other pulmonary and even non-pulmonary disorders.

The device is not intended for life support nor does it provide any patient monitoring capabilities.

1.2. Intended user

The Portable Mesh Nebulizer is intended for use by trained professionals and lay persons.

- The user should also be capable of understanding the general operation of this device and the content of this instruction manual.

1.3. Intended patient population

It is suitable for adults and children (≥ 6 months old) who require nebulizer therapy through nebulization mask or mouthpiece.

1.4. Intended environment

This Portable Mesh Nebulizer is intended for use in a medical facility, such as hospital, clinic and doctor's office, and in a general household environment.

1.5. Intended medical condition

Obstructive airway diseases (OADs), such as asthma, chronic obstructive pulmonary disease (COPD), and other pulmonary and even non-pulmonary disorders.

1.6. Intended clinical benefits

The advantages associated with the use of portable mesh nebulizer are that it does not require patient coordination between inhalation and actuation, thus making it suitable also for pediatric and geriatric patients.

2. Contraindications

- (1) Patient is allergic to aerosolized drugs.
- (2) Pentamidine powder.
- (3) Anesthetics agent.

3. Important Safety Notes


- Before use, ensure that there is no visible damage to the device or accessories. In case of any doubt, do not use the device and contact your retailer or the specified Customer Service address.
- Do not use health products or medicines containing essential oils for nebulization.
- You should always follow the instructions of your doctor regarding the type of medication to use, the dosage, and the frequency and duration of inhalation. Only use medication prescribed or recommended by your doctor or pharmacist.
- The use of this product for children and persons with special needs must be carried out under correct guidance and supervision.
- This unit is only used for specified purposes, only for nebulization. Do not use the device for any other purpose.
- Clean and disinfect the medication cup and accessories before using and before longer time storage.

- Please stop using the device if the components are damaged or the main unit is fully immersed in water accidentally.
- Keep the device away from your eyes when it is in use, as the nebulize medication could be harmful.
- Keep packaging material away from children (risk of suffocation).
- Do not use any additional parts that are not recommended by the manufacturer.
- If any serious incident has occurred in relation to the device, please report to the manufacturer and the competent authority of your member state immediately.

4. Working Principle

The working principle of the nebulizer is creating micro particles from a microporous mesh vibrating at ultra-sonic frequencies. A ceramic piezo-electric element is driving the ultra-sonic vibrations, powered by an electronic oscillation circuit. The micro particles are expelled from the microporous mesh and form an aerosol which can be directly inhaled. During inhalation, the particles are deposited in the patient's trachea, bronchus, alveoli and other areas. The medication is then absorbed through the mucosal tissue establishing the intended treatment.

5. Package Contents and Overview

Package Contents

A	Medication Cup	Including vibrating mesh
B	Main Unit	—
C	Nebulization Mask	Adult Mask
D	Nebulization Mask	Child Mask
E	Mouthpiece	—
F	USB Type-C Cable	—
G	Soft pouch	—
—	Instruction Manual	—

1	Medication Cup Cover
2	Spray Nozzle
3	Connection Electrode
4	Stand-by (Power button)
5	Type-C Connector
6	Mask Elastic Band
7	Drain Hole

6. Product Technical Parameters

Product Description	Portable Mesh Nebulizer
Model	NE-U300 REF NE-U300-EO
Power source	DC 2.4V, rechargeable Li-ion batteries
Charge	DC 5V, 1AAC adapter that complies with IEC 60601-1 It can be used for approximately 50 minutes on a single full charge.
Power consumption	<4.0 W
Nebulization rate	0.15 mL/min-0.90 mL/min
Aerosol output	Approx. 610 µg (2 mL, albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution)
Aerosol output rate	Approx. 101 µg/min (2 mL, albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution)
Percentage of fill volume emitted per minute	5.05 %
MMAD	4.12 µm Note: The MMAD test values of each device can vary, and the final product inspection report shall prevail.
GSD	1.72
Qr (<2 µm)	Approx. 7.4 %
Qr (2~5 µm)	Approx. 54.7 %
Qr (> 5 µm)	Approx. 37.9 %

Residual volume	≤0.3 ml (Measured in a vertical position)
The maximum temperature above ambient reached in the medication cup	≤30°C
Medication cup capacity	10 mL (max)
Appropriate medication quantity	2 ml minimum~10 ml maximum
Vibration frequency	130 kHz, Deviation ±10 %
Product dimensions/ Weight	Approx. 40 mm (D) x 64 mm (W) x 100 mm (H) / Approx. 117 g
Working environment	Temperature: +5°C ~ +40°C Relative Humidity: 15 % ~ 80 % R.H. Non-condensing Atmospheric Pressure: 860 hPa ~ 1060 hPa
Storage/delivery environment	Temperature: -20°C ~ +55°C Relative Humidity: 10 % ~ 93 % R.H. Non-condensing Atmospheric Pressure: 700 hPa ~ 1060 hPa
Mode of operation	Continuous operation
Classification of installation and use	Hand-held
Applied part classified	The classified of applied part is BF.

The maximum A-weighted sound pressure level	≤40 dB (A-weighted)
Service life *Conditions: 2 times per day	Main Unit: 5 years Medication Cup: 1 year Nebulization Mask: 1 year Mouthpiece: 1 year
Material	Enclosure: ABS Medication Cup: Inner layer: PC, outer layer: ABS Stand-by (Power button): Silicone Nebulization Mask: PVC Mouthpiece: PP
Applied parts	Stand-by (Power button), Nebulization Mask, Mouthpiece

Note:

- 1) According to the type of the classification of protection against electric shock: Internally powered ME equipment.
- 2) Enclosures classified according to degree of protection against the ingress of water and particulate matter as per IEC 60529: IP24 (2: Protected against solid foreign objects of 12.5 mm Φ and greater / 4: Protected against splashing water)
- 3) The equipment is not intended for use in an oxygen-rich environment and presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide.
- 4) The distribution of particle size is measured with Albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution. The particle size distribution and aerosol output may vary by a combination of product, medication and ambient conditions such as temperature, humidity and atmospheric pressure.
- 5) Performance may vary with the mesh and medication used such as suspensions or high viscosity. See the medication supplier's data sheet for further details.
- 6) The disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are

likely to be different from those stated for pediatric or infant populations.

- 7) Measured values reflect internal data obtained by NEXT GENERATION IMPACTOR (NGI), according to EN ISO 27427:2023.
- 8) The Nebulization rate is measured with albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution at $(23\pm 2)^{\circ}\text{C}$ with 2 ml of medication. It can vary with medication and ambient conditions.
- 9) Regarding the product technical parameters, unless otherwise specified, the following conditions are applied.
2 times a day for 10 minutes each time
Operating environment: $23\pm 5^{\circ}\text{C}$, 30 to 70 % RH, standard atmosphere (1013.25 hPa)
- 10) Higher than 70 % of the Initial Capacities of the Cells Carry out >2000 cycle
Charge: 1.0 C to 2.8 V Discharge: 1.0 C to 1.6 V Temperature: $25\pm 3^{\circ}\text{C}$.


★ Performance disclosures

The particle size and the aerosol output, aerosol output rate is measured according to the provisions of the EN ISO 27427:2023, Annex C (using breathing simulator) and D (using the multistage cascade impactor) separately.

The nebulization rate is measured with albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution at a temperature of $23\pm 2^{\circ}\text{C}$ and a relative humidity of 45-75 %. These values can vary with medication and ambient conditions.

The test solution is albuterol 0.1 % (m/m) concentration in 0.9 % sodium chloride solution. Test conditions of a temperature is $23\pm 2^{\circ}\text{C}$ and a relative humidity is 45-75 %. and the plot of cumulative size distribution of the results is as follows.

Aerodynamic Particle Size Distribution

Note:

- 1) The particle size distribution and aerosol output may vary by combination of product, medication and ambient conditions such as temperature, humidity and atmospheric pressure.
- 2) Using a solution, suspension, or emulsion different from that recommended by the manufacturer, in particular, a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer. See the medication supplier's data sheet for further details.
- 3) The effect of the following disclosures for nebulizer performance is based upon testing that utilizes adult ventilatory patterns and is likely to be different from those stated for pediatric or infant populations.
- 4) Essential performance: nebulization rate.
- 5) The proportion of aerosol particles with a diameter of 1 μm to 5 μm generated by the nebulizer is >50 %.

7. Installation and User Instructions


7.1. Power Supply

Use the USB Type-C cable (included in the package) and AC adapter (output: DC 5V, 1A complies with IEC 60601-1; not included in the package) for charging. It is recommended to use the optional AC Adapter Set Model NEB-ASET-3A or NEB-ASET-3AUK. (Please refer to "15. Configuration List".)

It takes approximately 2 hours to fully charge.

Charge the device as illustrated in figure below.

When the battery indicator light changes from flashes slowly to solid on, it's fully charged. (Please refer to "7.3. Indicator light".) Please disconnect the AC adapter.

Note:

If the battery indicator light is off, it may indicate that the AC adapter is incompatible or damaged or the equipment may have been damaged. In that case, please disconnect the AC adapter.

Warning:

- 1) At the end of expected service life, please dispose of this device and accessories according to the local environmental regulation, do not dispose together with the domestic refuse to avoid environment pollution.
- 2) Use only AC adapters with output DC 5.0V, 1A that comply with IEC 60601-1 to avoid electrical hazards.
- 3) Do not modify this equipment without authorization from the manufacturer.
- 4) Do not disassemble or repair this device, and do not disassemble or replace the battery of the manufacturer.


Caution:

- 1) When the battery capacity is insufficient, please power it using the USB Type-C cable.
- 2) Keep charging the device at least once per month during the storage period when this exceeds one month.
- 3) The device cannot operate while charging.
- 4) The main unit is powered by secondary lithium batteries which comply with IEC 62133-2.

7.2. Use Instructions Step

How to Use

- 1) Turn the medication cup cover counterclockwise and remove it.
- 2) Fill the medication cup with a minimum of 2 mL and a maximum of 10 mL of prescribed medication.
- 3) Install the medication cup cover.


Note:

- Before use, clean and disinfect the main unit, medication cup, and accessories. Refer to "8. Cleaning and Disinfection".
- Follow your doctor's advice/prescription to fill medication.
- Be careful not to overflow medication.


Caution:

- Fill the medication cup before turning the power on.
- The amount of medication should not be less than 2 mL and should not be more than 10 mL.
- DO NOT open the medication cup casually after filling the medication to avoid spilling.

4) Attach the Nebulization mask or mouthpiece to the medication cup.


a. Attach Nebulization mask


b. Attach mouthpiece

5) Place the Nebulization mask over mouth and nose or place the mouthpiece in your mouth.

a. Nebulization mask inhalation method

b. Mouthpiece inhalation method

6) Press the Stand-by (Power button) to start the nebulization. While using the device, sit upright and breathe calmly. Hold your breath for 2-3 seconds after inhalation to improve medication absorption.

⚠ Caution:

If too much force is applied to the nebulization mask or mouthpiece, the medication cup may come off from the main unit.

During use, hold the main unit and the medication cup firmly and be careful not to separate them.


7) After the inhalation, press the Stand-by (Power button) to stop the nebulization. The nebulization automatically stops after 30 minutes.

8) After use, remove the medication cup cover and discard the medication.

a. Turn the medication cup cover counterclockwise, then pull it upward to remove.


b. Pour out any residual medication by tilting downward.

9) Clean and disinfect the components. Refer to "8. Cleaning and Disinfection".

Caution:

- 1) When installing the medication cup, do not tilt, make sure the whole medication cup is connected securely.
- 2) Keep the device vertical during the operation, no shaking during nebulization.
- 3) Holding your breath for a few moments before starting the treatment can enhance the effectiveness of nebulization therapy. Sit down and keep calm and relax, do not inhale too fast during the treatment process.
- 4) Replace the nebulization mask, or mouthpiece if the inhalation method is changed. Please clean the medication cup after use (Refer to "8. Cleaning and Disinfection").
- 5) If you notice that liquid has built up around the nebulization nozzle and vibrating mesh after using the device, it may affect how well the medication is nebulized. In this case, please stop the nebulization process. Next, carefully remove the inhalation mask and any other accessories. Use a clean medical gauze to gently wipe away any leftover liquid. Please be careful not to touch the central spray area of the vibrating mesh, as this could cause damage.
- 6) For questions on the process of nebulizing, please contact manufacturer or Distributor.

7.3. Indicator light

Operation indicator solid on	Operating
Operation indicator flashes rapidly 10 times	Low medication, and the device will automatically shut down
Operation indicator flashes slowly 10 times	30-minute nebulization timer ends, and the device will automatically shut down
Battery indicator flashes rapidly 10 times	Low battery
	Battery depleted, and the device will automatically shut down
Battery indicator flashes slowly	Charging
Battery indicator solid on	Fully charged

Caution:

(1) Definition of “flashing rapidly”: Light on for 0.5 seconds, off for 0.5 seconds.

Definition of “flashing slowly”: Bright → dim for 1 second, dim → bright for 1 second (alternating without turning off)

(2) When a low medication warning appears and the operation indicator flashes rapidly 10 times

Adding the medication will resume nebulization;

Pressing the Stand-by (Power button) alone will not forcibly start the nebulization.

EN

7.4. Software version

WHQW8 V1.0

8. Cleaning and Disinfection

Parts which can be cleaned and disinfected:	Medication cup, Nebulization masks and Mouthpiece
Advice:	Reprocessing procedures have only limited implications to this product. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. From the processing side there is no maximum number of allowable reprocessing. The device should not be reused if there are signs of material degradation, e.g. cracks on the product, brittle, other change in the material, etc.
Cleaning and Disinfection Instructions	
Preparation at the Point of Use:	Remove visible residues on the device with cold water (<40°C) immediately after use, if applicable. Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residues which may influence the result of the reprocessing process.
Transportation:	Safely store the device in a humid environment and transport it to the reprocessing area to avoid any damage and contamination to the environment.
Preparation for Decontamination:	The devices must be reprocessed in a disassembled state, as far as possible.
Pre-Cleaning:	Perform manual pre-cleaning, until the devices are visually clean. Submerge the components in a cleaning solution. Clean the surfaces with a soft lint-free cloth.

Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better reproducibility and standardization, and for personnel protection.

Automated Cleaning:

Use a washer-disinfector meeting the requirements of the ISO 15883 series.

The instruments are placed on the cleaning rack of the WD. The instruments in the washer-disinfector are arranged in such a way that there is no rinsing shadow and the water drains off quickly.

Start the program:

Cleaning:

- 4 min. pre-washing with cold water (<40°C)
- emptying
- 5 min. washing with a mild alkaline cleaner at 55°C
- emptying
- 3 min. rinsing with warm water (40°C)
- emptying
- 5 min. intermediate rinsing with warm water (40°C)
- emptying

The automated cleaning processes have been validated by using 0.5 % neodisher MediClean forte.

Note fresh de-ionized water should be used for the automated cleaning process. Water quality refers to EN 285.

Manual Cleaning:

Manual cleaning of the parts can be performed under the following conditions of: 0.5 % neodisher MediClean forte, 10 minutes immersion time at 30°C.

If other detergents are used, select immersion times, concentrations and temperatures of the solutions as instructed by the manufacturer of the detergent.

In that case, it is recommended to use a detergent containing alkalinity donors, enzymes, and surfactants.

- a. Wear disposable protective gloves.
- b. Put the parts in a suitable cleaning basin.
- c. Create a detergent solution as instructed by the manufacturer of the detergent.
- d. Completely immerse the parts in the detergent solution. Make sure that there are no air bubbles on the products while immersed.
- e. Dip a clean soft lint-free cloth in the detergent solution, and gently wipe all surfaces and lumen of the parts until these are visually clean. Pay special attention to the internal parts. When the cloth becomes obviously contaminated, discard it, and use a new cloth as needed. Then flush the lumen of the parts with the detergent solution by using a syringe. Repeat the step for 5 times. The cleaning process takes about 10 minutes.
- f. Remove the parts from the cleaning basin. Rinse the device under running water for 30 seconds at room temperature.
- g. Put the cleaned and rinsed parts into a suitable clean basin with fresh de-ionized water. Thoroughly rinse and flush the lumen of the parts with the water for 5 minutes.
- h. Dry the parts thoroughly with a lint free towel.

Automated Disinfection Program (Glutaraldehyde Disinfectant):

Start the program:

- 5 min. disinfecting with Belimed Protect™ Glutaraldehyde Disinfectant at 55°C
- 1 min. neutralizing with water at ambient temperature
- emptying
- 1 min. neutralizing with water at ambient temperature
- emptying
- 8 min. drying at 50°C

Manual Disinfection:

Following parameters are used:

Disinfectant: MetriCide OPA Plus (* active ingredient: ortho-phthalaldehyde)

Concentration: ready-to-use solution, undiluted

Temperature: ambient temperature

Soaking time: ≥5 min.

- a. Immerse the parts in the disinfectant solution.
- b. Wipe the surface of the parts with a sterile soft lint-free cloth to remove any air bubbles on the surface.
- c. After completion of soaking time, remove the parts from the disinfection solution.
- d. Rinsing the parts after manual disinfection with sterile water.
- e. Wipe and dry the parts surfaces with a sterile lint-free cloth until the surface is free of water stains.

The disinfection conditions described above are only recommended parameters, please refer to instructions for use of the disinfectant agent for actual use parameters and cautions.

Disinfection:

Drying:	Drying the parts thorough drying cycle of washer/disinfector. If needed, (additional) manual drying can be performed through lint free towel. Blow-dry cavities of products by using sterile compressed air. For manual reprocessing, drying takes place following the rinsing step after manual disinfection. Shake the medication cup well to remove any liquid that has entered the cup through the drain hole at the bottom of the medication cup. Wipe and dry the surface with a sterile lint-free cloth until the surface is free of water stains. Blow-dry cavities of parts by using sterile compressed air.
Functional testing, Maintenance:	Visual inspection for cleanliness of the parts. All parts should be checked again for dryness. After cleaning and disinfection, a thorough inspection and maintenance ensures that the parts are fit for use. - Check that the parts have no dents, cracks, deformations, scratches, etc. - Check all markings on the product for clear visibility. Discard and replace any components when necessary. Do not use the parts with following defects: material deformation, cracks on the product, parts which have become brittle or other changes in the material, etc. Test and maintain the parts according to manufacturer's instruction.
Functional test:	1. Visually inspect each part of the system for cracks or damage and replace if any defects are visible. 2. Pour 1-10 mL of normal saline (0.9 %) into the nebulizer. 3. Press the Stand-by (Power button) and check if the operation indicator is on and aerosol is visible. 4. Press the Stand-by (Power button) again to turn the system off and check if the operation indicator is off. Discard any remaining liquid before patient use.

Storage:	Store the device and the processed parts in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.
Additional Instructions: None	
It is the responsibility of the user to ensure that the reprocessing processes including resources, materials and personnel are adequate to reach the required results. State-of-the-art and often national laws require that these processes and included resources be validated and maintained properly.	

**Warning:**

- 1) Please don't use disinfectants containing Benzalkonium bromide or house bleach.
- 2) Be careful not to damage the vibrating mesh when disinfecting and cleaning. Do not use medical gauze to wipe the residues in the vibrating mesh, let the parts air dry.
- 3) Never use Sodium hypochlorite, Hypochlorous acid or Quaternary ammonium compound to disinfect.

Note:

- 1) The equipment's water resistance classification is IP24. This main unit cannot be washed by flowing water prevent the ingress of water at any time.
- 2) The disinfectant residues should be cleaned thoroughly to assure safe use next time.
- 3) Before cleaning and disinfection shall be carried out the power source must be disconnected.
- 4) After cleaning and disinfection reassemble the device according to 7.2.

9. Storage and Maintenance

- 1) The storage and transportation conditions of the products are detailed in "6. Product Technical Parameters".
- 2) Keep the device out of reach of unsupervised infants and children. The small parts of this device may be swallowed by the infants and children.
- 3) Prevent pets and pests from damaging this device.
- 4) Dry the parts immediately after cleaning and disinfection. Store the device and the components in the environment that meets the requirements, be careful to avoid mechanical shocks like dropping the device.
- 5) Screw off the cap and discard the residual liquid from the medication cup. Do not reuse.
- 6) Please wipe the outside of the Main Unit by a soft cloth.
- 7) The device should be stored in a place with no corrosive gas and good ventilation, please avoid severe shocks during transportation.
- 8) Do not use benzene, diluent and flammable chemicals to clean the product.

10. Trouble Shooting

Item	Trouble	Possible cause/solution
1	Does not work when switched on.	• If the battery indicator flashes or no indicator flashes or is lit, recharge the device. Refer to "7.1. Power Supply". • If the operation indicator flashes rapidly, the medication level is low. Fill the medication cup and try again. • Clear the clogged medication in the vibrating mesh, then try again. • If the device does not nebulize after checking the above causes, contact your OMRON retail outlet or distributor.
2	Low nebulization rate.	• Check if the medication cup is filled with right medication, which should be water-solubility, non-corrosive medication. • Check if the medication cup is filled with appropriate volume. • Tilt Main Unit, so that the medication can contact with the vibrating mesh. • Re-assemble the medication cup correctly and restart the device. • Clean the medication cup and vibrating mesh. If it still cannot be used after cleaning, please check whether the vibrating mesh is broken or not. • Clear the clogged medication in the vibrating mesh and restart the power.
3	The device turns off immediately after switching on	• Install the medication cup correctly, then try again. • Tilt the device slightly forward to bring the medication into contact with the vibrating mesh. • Clear the clogged medication in the vibrating mesh, then try again. • If the device does not nebulize after checking the above causes, contact your OMRON retail outlet or distributor.

4	The mesh is clogged with medication.	• Take off the medication cup and the medication cup cover. Wash with warm water and air dry.
5	The medication cup assembly leakages.	• Re-assemble the medication cup correctly. • Purchase the new medication cup. Contact your OMRON retail outlet or distributor.
6	Residual medicine liquid in the medication cup.	• Always clean the medication cup after using. Follow the cleaning and disinfection procedure.
7	Can not charge.	• Make sure the USB Type-C cable and the AC adapter are correctly connected. • Purchase the new USB Type-C cable or AC adapter. Contact your OMRON retail outlet or distributor.

11. Note & Warning

11.1. Note and suggestion

1. This device is a medical device, please read the user manual before using it.
2. Please use original accessories, warranty service is not provided for damage caused by accessory beyond our list.
3. Use of a non-conforming adapter may cause equipment damage.
4. Please refer to the user manual when there is a problem or contact the after-sales service unit for maintenance.
5. Please clean and disinfect the medication cup before using it for the first time. You can refer to "8. Cleaning and Disinfection".
6. This device is for medication nebulization, not for humidification, please avoid the use of distilled water, it may cause lower nebulization rate.
7. Please clean the medication cup empty before storing it.
8. Please make sure all the accessories are assembled correctly before using.
9. Please use the accessory individually to avoid cross infection.
10. Please keep the device in vertical position when nebulizing.
11. Don't use the unit under inflammable gas environment or near a heating device or near an open flame.
12. Don't use the unit near high frequency products or electronic products.
13. Do not use a microwave oven, oven, blow dryer or other house applications to dry nebulizer and accessories.

14. Technical description is included in this user manual.
15. The nebulizer can only be used by operators who can understand this user manual.
Children should use under adult supervision only.
16. When the product is taken out from -20°C, it should acclimatize to room temperature for 2 hours before use.
17. When the product is taken out from +55°C, it should acclimatize to room temperature for 2 hours before use.
18. The charging interface is only used for power supply and cannot be used to connect to other devices.
19. The patient can be the intended operator and all functions of the devices can be used by patients.
20. When the patient is an intended operator, all the accessories of the nebulizer cannot be serviced or maintained while it is in use.
21. Maintenance for the device by the patient is limited to cleaning and disinfection. Please refer to "8. Cleaning and Disinfection".
22. Performance information provided in this instruction manual is in accordance with EN ISO 27427:2023. It may not apply to drugs supplied in high viscosity form. In such cases, please consult your doctor regarding the type of drug you are taking.
23. The highest rated operating altitude:1300m (which corresponds to an ambient pressure of 872 hPa)
24. The nebulizer device and the accessories (nebulization mask and mouthpiece) comply with the requirements of ISO 18562-2, ISO 18562-3, ISO 18562-4.
25. Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

- | |
|---|
| 26. The materials used in the components could not be compatible with solutions/suspensions/emulsions that have not been evaluated.
The medicines that have been evaluated are: Ambroxol hydrochloride solution for inhalation and Salbutamol sulfate solution for inhalation. |
| 27. Check the power-on status after reassembling the device. |
| 28. Please only use or purchase original parts or accessories. |
| 29. These specifications are subject to change without notice. |

11.2. Warning

1. Please stop using the device if you feel uncomfortable and please consult your doctor.
2. Volatile oils are not allowed, this may cause damage to mesh.
3. Water-soluble medication and saline dilution medication are allowed but may cause Bronchospasm.
4. Please do not use medicines containing esters, oils or suspended particles, including herbal extracts. It is recommended to use the device according to the doctor's instructions.
5. Do not service and maintain when the nebulizer is in use.
6. Don't modify the equipment without authorization of the manufacturer. Otherwise, it may cause damage to the Main Unit or be harmful to the user or patients.
7. Replacement of lithium batteries when incorrectly replaced would result in an unacceptable risk.
8. This device is intended for single patient use and can be reused by a single patient.
9. Do not use mobile (cellular) telephones and other devices(such as MRI, diathermy, electrocautery, RFID and electromagnetic security systems) which generate strong electrical or electromagnetic fields, near the medical device.
10. The manufacturer or manufacturer's representative can assist the lay operator or lay responsible organization to set up, use or maintain this device when needed. In case of any unexpected operation or events, contact your OMRON retail outlet or distributor.

11. Battery disposal

The empty rechargeable batteries must be disposed of through dedicated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries according to local environmental regulations.

General disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its life. Dispose of the device at a suitable local collection or recycling point in your country in accordance with EC Directive -WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

The NE-U300 rechargeable battery is not intended to be disassembled and disposed of by users.

12. The USB Type-C cable should be kept away from children or out of the reach of children to avoid strangulation.

13. Consult your doctor for help or contact the manufacturer if any changes in the performance of this device.

14. Be careful with the small parts of this device because of the risk these may be swallowed by the infants or children.

15. Applied parts are not intended to supply heat to a patient.

16. The nebulizing system is not suitable for use in an anesthetic breathing system or a ventilator breathing system.

17. The types of liquid used in the device are designed for nebulization.

18. The accessible materials used in the nebulizer are safe for regular people. For very few operators with extreme skin sensitivity, if any skin discomfort occurs during the use of the nebulizer, please immediately stop using and seek medical advice.

19. Do not use the nebulization mask or mouthpiece or medication cup on multiple patients. (Single Patient Use Only). Re-use of single-patient-use parts should only be according to the advice from doctor.
20. The surface temperature of the applied parts exceeds 41°C when test in both normal condition and single fault condition. (Outside plastic enclosure near battery: 47.2°C. Power switch: 46.3°C. Medication cup: 45.0°C.) Contact time of applied parts during operation: 1 min. $\leq t < 10$ min.
21. Do not bring magnetic cards, magnetic media, precision instruments, or other devices that are susceptible to magnetism close to the main unit. (The built-in magnet may affect these items and cause damage to magnetically recorded data or lead to malfunctions.)

12. EMC declaration

1. This device complies with the requirement of electromagnetic compatibility in EN 60601-1-2:2015 + A1:2021.
2. The user needs to install and use the device according to electromagnetism compatibility information provided in this instruction manual.
3. Portable and mobile RF communication devices and some household appliances, such as mobile, interphone, microwave oven, dry blower, may influence Portable Mesh Nebulizer performance, so Portable Mesh Nebulizer should be kept away from them during use.
4. Guidance and manufacturer's declaration stated in the appendix.
5. Do not use this device in the vicinity of HF surgical equipment.

1. With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices.
 2. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation.
 3. Medical devices should also not interfere with other devices.
 4. In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN 60601-1-2:2015 + A1:2021 standard has been implemented. This standard defines the levels of immunity to electromagnetic interference as well as maximum levels of electromagnetic emissions for medical devices.
 5. Medical devices manufactured by FEELLIFE HEALTH INC. conform to this EN 60601-1-2:2015 + A1:2021 standard for both immunity and emissions.
- Nevertheless, special cautions need to be observed:

Note:

The EMISSIONS characteristics of this equipment make it suitable for use in a residential environment (for which CISPR 11 class B is required).

**Warning:**

1. The use of accessories not specified in this manual may result in increased emissions or decreased immunity of the device.
2. The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.
3. PORTABLE RF communications equipment(including peripherals as such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could occur.
4. Do not use mobile (cellular) telephones and other devices which generate strong electrical or electromagnetic fields, near the medical device. This may result in incorrect operation of the unit and create a potentially unsafe situation. Recommendation is to keep a minimum distance of 7 m. Verify correct operation of the device in case the distance is shorter.

Details of cable:

Name	Type	Length (m)	Cable Shielded
USB Type-C cable	USB Type-C	1.0	UnShielded

Guidance and Manufacturer's declaration – electromagnetic emissions		
This Portable Mesh Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user of this Portable Mesh Nebulizer should assure that it is used in such an environment.		
Immunity test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This Portable Mesh Nebulizer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This Portable Mesh Nebulizer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's declaration – electromagnetic immunity			
This Portable Mesh Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user of this Portable Mesh Nebulizer should assure that they are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line to line	± 0.5 kV, ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, Voltage interruptions IEC 61000-4-11.	0 % U_T , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T , 1 cycle and 70 % U_T , 25/30 cycles at 0° 0 % U_T , 250/300 cycle	0 % U_T , 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T , 1 cycle and 70 % U_T , 25/30 cycles at 0° 0 % U_T , 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Portable Mesh Nebulizer requires continued operation during power mains interruptions, it is recommended that the Portable Mesh Nebulizer be powered from an uninterruptible power supply or a battery.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz & 60 Hz	30 A/m 50 Hz & 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic fields IEC 61000-4-39	CW 8 A/m for 30 kHz Pulse modulation 2.1 kHz, 65 A/m for 134.2 kHz Pulse modulation 50 kHz, 7.5 A/m for 13.56 MHz	CW 8 A/m for 30 kHz Pulse modulation 2.1 kHz, 65 A/m for 134.2 kHz Pulse modulation 50 kHz, 7.5 A/m for 13.56 MHz	Power frequency magnetic fields / Proximity magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance & Declaration – Electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that they are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment –guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands and amateur radio bands	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands and amateur radio bands	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Portable Mesh Nebulizer, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 385 MHz to 5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014+A1:2020)	10 V/m 80 MHz to 2.7 GHz 385 MHz to 5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014+A1:2020)	

13. Signs and Symbols

General symbol	Description / title	General symbol	Description / title
	Date of manufacture		General warning sign
	Manufacturer information		Refer to instructions for use (Background: blue)
	Authorized representative in the European community		Caution
	Temperature limit		Serial Number
	Reference		Humidity limitation

	Medical device		Atmospheric pressure limitation
	Stand-by		When this device life expires, the end users should discard this device according to the requirements from the local environment protection authority.
	Direct Current	IP	Ingress protection degree provided by IEC 60529.
	Recycle mark X: Material number Y: Material abbreviation Refer to 97/129/EC for more information.	CE 0197	Notified Body for CE mark
	Recycling instruction for packaging elements	UDI	Unique device identifier
	Indoor use		Energy efficiency level 6

	Do not use if package is damaged		Double insulation
	Batteries should be collected separately		Batch code
	Need for the user to consult the instructions for use		Alternating current
	Type BF applied part		

14. Warranty

Thank you for buying an OMRON product. This product is constructed of high-quality materials and great care has been taken in its manufacturing. It is designed to give you a high level of comfort provided that it is properly operated and maintained as described in the instruction manual.

This product is guaranteed by OMRON for a period of 3 years after the date of purchase (except for replacement parts). The proper construction, workmanship and materials of this product is guaranteed by OMRON. During this period of guarantee OMRON will, without charge for labor or parts, repair or replace the defective product or any defective parts.

The guarantee covers only products purchased in Europe, Russia and other CIS countries, Middle East and Africa.

The guarantee does not cover any of the following:

- a. Transport costs and risks of transport.
- b. Costs for repairs and / or defects resulting from repairs done by unauthorized persons.
- c. Periodic check-ups and maintenance.
- d. Failure or wear of accessories (including Medication Cup) or other attachments other than the main device itself, unless explicitly guaranteed above.
- e. Costs arising due to non-acceptance of a claim (those will be charged for).
- f. Damages of any kind including personal caused accidentally or from misuse.

Should guarantee service be required please apply to the dealer whom the product was purchased from or an authorized OMRON distributor. For the address refer to the product packaging / instruction manual or to your specialized retailer. If you have difficulties in finding OMRON customer services, visit our website (www.omron-healthcare.com) for contact information.

Repair or replacement under the guarantee does not give rise to any extension or renewal of the guarantee period. The guarantee will be granted only if the complete product is returned together with the original invoice / cash ticket issued to the consumer by the retailer. OMRON reserves the right to refuse the guarantee service if any unclear information has been given.

15. Configuration List

Optional Medical Accessories

Product description	Model	REF	Remarks
Medication Cup	MAW8	NEB-MEC-3AE	—
Nebulization Mask	L2	NEB-MSLP-3AE	Adult Mask
	S2	NEB-MSMP-3AE	Child Mask
Mouthpiece	M1	NEB-MP-3AE	—

Other Optional Accessories

Product description	Model	Remarks
AC Adapter Set	NEB-ASET-3A	USB Type-C cable (UTC021) C type AC Adapter (SINGOF-5E-050100)
AC Adapter Set (UK)	NEB-ASET-3AUK	USB Type-C cable (UTC021) BF type AC Adapter (SINGOF-5K-050100)

16. Disclaimer Clause

Please read the user manual before using this device. The manufacturer will not take any responsibility in case of damage caused by improper use of this device.

Please use or purchase original parts or accessories. The manufacturer does not take responsibility for the buyer or third parties for any damage or loss intentionally or unintentionally caused by improper use.

On the request for warranty service, please present your warranty card filled with purchase date and seal (with the store and address). Any repair service out of the scope of warranty will be charged accordingly.

Manufacturer 	FEELLIFE HEALTH INC. Room 2501, 1903, 2002, Building A, No.9 Furong Road, Tantou Community, Songgang Subdistrict, Bao'an District, Shenzhen, 518104, Guangdong, P.R. China Tel: +86-755-66867080 E-mail: info@feellife.com Website: https://www.feellife.com
EU-representative EC REP	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595 AA, The Hague, The Netherlands Tel: +31644168999 Email: info@lotusnl.com
Distributor	OMRON HEALTHCARE EUROPE B.V. Wegalaan 73, 2132 JD Hoofddorp THE NETHERLANDS www.omron-healthcare.com
Importer in EU	
Production facility	FEELLIFE HEALTH INC. Room 2501, 1903, 2002, Building A, No.9 Furong Road, Tantou Community, Songgang Subdistrict, Bao'an District, Shenzhen, 518104, Guangdong, P.R. China Tel: +86-755-66867080 E-mail: info@feellife.com Website: https://www.feellife.com
UK Responsible Person	KINGSMEAD SERVICE LIMITED 19 MEZZANINE FLOOR 19-21 CRAWFORD STREET LONDON ENGLAND W1H 1PJ
Importer in the United Kingdom	OMRON HEALTHCARE UK LTD. SEEBECK House 1, SEEBECK Place Knowlhill, Milton Keynes, MK5 8FR, United Kingdom www.omron-healthcare.com/distributors
Subsidiaries	OMRON MEDIZINTECHNIK HANDELSGESELLSCHAFT mbH www.omron-healthcare.com/distributors
	OMRON SANTÉ FRANCE SAS www.omron-healthcare.com/distributors

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