



HYPERVOLT

Featuring Quiet Glide™ Technology

As far back as I can remember, I have always had a passion for sports, fitness and science. From my own experience as an athlete, I understood that maintaining a healthy body is essential to an athlete's performance. Several years ago, I set out to create a brand whose products would improve athletic performance by enhancing mobility, recovery, and preventing injury. In order to develop products of this caliber, I collaborated with some of the world's best athletes and trainers to develop products that would meet their standards. I would go on to name the brand after our first product – the HYPERICE: the most advanced portable ice compression device on the market. Since then, we have gone on to create a series of products that elite athletes from all over the world have adopted as part of their daily routine. Our goal is to take these innovative products and put them in the hands of anyone, regardless of fitness level, looking to improve their overall health.

We believe this vibration massage device will appeal to anybody suffering with pain, soreness, and stiffness. The Hypervolt is quickly becoming one of the most revered vibration massage devices today.

*AK* Anthony Katz, Founder

**THE HYPERVOLT**

# Power in the palm of your hand.

The Hypervolt is a cordless state-of-the-art vibration massage device that relaxes sore stiff muscles, increases circulation and range of motion, and improves the overall health of the body's soft tissue.



**CORDLESS / RECHARGEABLE**

24 V Lithium-Ion Battery



**LIGHTWEIGHT DESIGN**

2.5 lbs / 1.1 kg

## SPECIFICATIONS



- 1 Powerful Brushless High-Torque Motor with Quiet Glide™ Technology
- 2 Removable Battery Release
- 3 Rechargeable Lithium-Ion Battery (2+ Hrs of Use Per Charge)

*\*Additional batteries available at hyperice.com*

- 4 3 Speed Settings for Optimal Pressure
- 5 Cordless and Lightweight (2.5 Lbs)

*\*Includes AC wall charger*

- 6 ON / OFF Power Switch
- 7 Interchangeable Applicators (x4)

### Vibration Frequency:

Level 1 – 30Hz, Level 2 – 40Hz, Level 3 – 53Hz

### Weight:

2.5 lbs / 1.1 kg

### Charge:

110-240V 50/60Hz – .5A – 24VDC-1A

## USAGE INSTRUCTIONS

PLEASE READ PRIOR TO USE

### CHARGING

- Fully charge battery for up to six hours before first use either attached or detached from the device.
- To charge, connect the DC end of the supplied 24V adapter to the battery's charging port and plug AC side cable into wall outlet.
- LED band on battery lights pulsate to show battery level and indicate active charging.
- With battery attached, LED lights on unit flash to show battery level and indicate active charging.
- Full charge is indicated when all LED lights stop flashing and remain fully illuminated.
- Battery may be recharged at any time and at any battery level.
- It is not recommended to fully deplete the battery to the red LED level.
- The average run-time is 2+ hours depending on speed level and pressure applied during use.

### OPERATING

- Align and attach battery into base of unit – ensure it is securely locked in place.
- Power device on by moving the power switch on the bottom of the battery to the ON position. Power device off by moving the power switch to the OFF position.
- With power switch in ON position, press the power button once to turn device on to level one, twice to turn device to level two and three times for level three. Pressing the power button again will turn the device off. The blue LED lights will indicate the corresponding power level.
- With unit off, fully insert applicator by pressing in. Remove by pulling it out.
- Massage desired body part at the desired speed applying pressure as needed.
- To remove the battery, lightly press the battery release button with thumb. Firmly pull downwards on the battery / handle with the other hand.

### MAINTENANCE, CLEANING & STORAGE

- Use a slightly damp towel to clean device and a dry with a soft cloth.
- For storage or travel, remove battery and store in supplied carton box.

## DISCLAIMER

## DISCLAIMER

### WARNINGS & CAUTION

- Do not immerse in water. Keep away from liquids or heat sources. Keep ventilation ports free from dust and debris.
- Do not remove screws or attempt to disassemble.
- Do not operate continuously for more than one hour. Allow device to rest for 30 minutes before re-using.
- Unplug the charger from device battery after charging or prior to use.
- For adult use only. Do not use if injured. Consult your doctor before using this product.

### SAFETY INSTRUCTIONS

**DANGER: TO REDUCE RISKS OF ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:**

- For adult use only.
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds per region.
- Use the Hypervolt only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body.
- Use only the applicators that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort.
- Keep fingers, hair or any other body part away from the shaft and backside of the applicator as pinching may occur.
- Do not place any objects into the ventilation holes of the Hypervolt.
- Do not immerse in water or allow water to enter the ventilation holes of the Hypervolt.
- Do not drop or misuse the Hypervolt .
- Only recharge with supplied 24VDC HYPERICE charger.
- Carefully examine device and battery before each use.
- Do not tamper or alter Hypervolt in any way.
- Never leave the Hypervolt operating or charging unattended.



**PLEASE DO NOT USE THE HYPERVOLT, OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR IF ANY OF THE FOLLOWING APPLY:**

Pregnancy, diabetes with complications such as neuropathy or retinal damage, wear of pace-makers, recent surgery, epilepsy or migraines, herniated disks, spondylolisthesis, spondylol-ysis, or spondylosis, recent joint replacements or IUD's, metal pins or plates or any concerns about your physical health. Frail individuals and children should be accompanied by an adult when using any vibration device. These contra indications do not mean that you are not able to use a vibration or exercise device but we advise you to consult a doctor first.

Ongoing research is being done to the effects of vibration exercise for specific medical disorders. This will likely result in the shortening of the list of contra indications as shown above. Practical experience has shown that integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advice of and in accompaniment with a doctor, specialist or physiotherapist.

### WARRANTY

**IF WITHIN ONE YEAR FROM THE DATE OF PURCHASE, THIS PRODUCT FAILS DUE TO A DEFECT IN MATERIAL OR WORKMANSHIP, HYPERICE, INC. WILL REPAIR OR REPLACE THE PRODUCT, OR NECESSARY COMPONENTS, FREE OF CHARGE.**

**THIS WARRANTY EXCLUDES:**

- (a) damage caused by accident, abuse, mishandling, or transport;
- (b) units subjected to unauthorized repair;
- (c) units not used in accordance with Hyperice care instructions;
- (d) damage exceeding the cost of the product;
- (e) deterioration of the delivered product resulting from abnormal storage and/or safeguarding conditions on the client's premises, and
- (f) failure to provide the dated proof of purchase.

Some states and countries do not allow a limitation of damages, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, and but other rights may vary from country to country and from state to state. This warranty service is available through the operating policies and procedures of Hyperice Inc. and all Hyperice international distributors. If you are a client, please contact the appropriate Hyperice business entity, dealer or retailer, from whom you directly purchased the product for warranty and return authorization procedures. This warranty is initiated and executed by Hyperice Inc.

## OTHER PRODUCTS

# Vyper™ 2.0

## High-Intensity Vibrating Fitness Roller

The Vyper 2.0 is the next generation of the world's most powerful vibrating fitness roller. The award winning Vyper 2.0 features 3 speeds of high-intensity vibration and a dual-zone smooth and grooved German engineered exterior. The Vyper 2.0 allows you to warm up, activate and recover faster and more effectively than any other roller.

### Available Colors

Black, Electric Green, Grey Camo



## OTHER PRODUCTS

# Hypersphere™

## High-Intensity Vibrating Fitness Ball

Compared to other SMR / massage balls, the Hypersphere helps to pinpoint and release trigger points faster, deeper, and less painfully using high-intensity vibration. The Hypersphere is portable and only 5" in diameter, which helps to release tension in targeted areas better than any roller.

### Available Colors

Electric Green, Stealth Black



## OTHER PRODUCTS

# Venom™

## Wearable Vibration + Heat

The Venom is a cutting-edge wearable device that uses heat + vibration to warm up, loosen, and relax sore/stiff muscles. Manufactured using the industry's best technology, the digital touchscreen allows you (the user) to customize the temperature, vibration pattern, and timer for the optimal massage.

### Available Devices

Back, Leg, Right Shoulder, Left Shoulder



## OTHER PRODUCTS

# Hyperice® ICT

## Ice Compression Technology w/ Patented Air-Release Valve

The ICT is an easy-to-use, portable 2-part cryotherapy device. Combining a Plush™ neoprene Compression Sleeve + Ice Cell with a patented air-release valve, the ICT molds the ice to your designated body part perfectly for an optimal "ice cast" treatment.

### Available Devices

Back, Left / Right Shoulder, Knee, Utility (Ankle, Elbow, Shin, Wrist)



## OTHER PRODUCTS

# Relaxroll®

## Superior Fitness Rollers

Relaxroll is a line of eco-friendly rollers uniquely designed to fit every contour of the body for a more comfortable roll.

**Relaxroll** - Contours to the curve of the back and lower extremities

**AchillX** - Ergonomically shaped for the calf, foot, neck, forearm, and Achilles

**SphereX** - Pin-point areas in the back, hips, glutes, shoulders, and hamstrings



## OTHER PRODUCTS

# Hypersoothe™

## Essential Oil-Infused Soothing Cream for Muscles and Joints

Hypersoothe is an all-natural, lightweight, topical cream infused with essential oils and natural plant extracts to soothe the body's muscles and joints. Apply this fresh fragrant topical before/after physical activity, throughout the day, or during a massage treatment.

Peppermint Oil to cool, Aloe and Nettle to soothe, Lavender and Chamomile to calm, Coconut Oil to moisturize and nourish, Sweet Almond Oil to soften skin, and Tea Tree Oil to disinfect and purify.





15440 Laguna Canyon Rd, Suite 230  
Irvine, California 92618  
(949) 565-4994  
[hyperice.com](http://hyperice.com)